

Fudo Mitten Side



Features a side zipper

- Opens wide from below

Wrist → Fingertip



Easy to insert hand

Protective Cover

- Opens from the top as well

Fingertip → Wrist



Finger Checking

Flap

Ambidextrous

Blue

Pink

Green

Can be opened and closed from the top and bottom



Side zipper with a covered flap and a hidden "Idea Hook" closure, designed to conceal the zipper tab.



Side Zipper

Palm Side

Back of the Hand


Features



Adjustable Belt



Designed with a cover over the belt to make it difficult to remove.



Replaceable "Idea Hook" belt cover



[Optional Replaceable Cover]

Standard "Idea Hook" type Belt Cover

Blue

Pink

Green

Material

Glove Fabric : Cotton 100%

Inner Fabric : Polyester 65% · Cotton 35%

Mesh · Trim Belt : Polyester 100%

Washing Instructions

Make sure hooks and buckles are closed before washing.
Washing them open or detached can lead to wear or damage.



Size

	M	L
Wrist	12~18cm	15~21cm
Hand Length	14~18cm	18~22cm



This product is used temporarily to support medical and caregiving procedures when there is a risk of the patient engaging in dangerous behaviors, such as self-harm or pulling out catheters, or when they try to tamper with their diaper.

"Idea Hook"



[How to unhook]
Support the back of the hook with your finger, and press the center with your fingernail.



[How to close]
Place both hands on the edges of the hook over the flap and press firmly until you hear a 'click'.



<https://www.taketora-web.com/>

Taketora Co., Ltd.

[Human Care Department]

8F Yokohama Media Business Center, 2-23 Otamachi,
Naka-ku, Yokohama-shi, Kanagawa 231-0011, Japan
TEL +81-45-680-0390 FAX +81-45-680-0399